

Packaged Fumigation Transfer Hatch.

BSL3/BSL4 Containment Solution



- Choice of fumigants HPV or Formaldehyde
- PLC control with colour touch screens
- Onboard HEPA filters and pressure relief filter
- 316 Stainless construction
- APR doors (mechanical or inflatable seal options)



PBSC Bio decontamination chambers are in use throughout the world – Years of expertise have been distilled down into a cost effective standard package. Easy to install, easy to validate

For more information contact:

PBSC Ltd. Bradley Junction Industrial Estate – Leeds Road – Huddersfield UK

Tel +44 (0)1484 354500 Email info@pbsc.co.uk Web www.pbsc.co.uk

